

DI 18 basic

Dispersing device

- An affordable low cost high value unit for simple mixing tasks
- The wide speed range makes it possible to work at high circumferential speed
- Accepts plug-in dispersing tools

Volume range (H ₂ O)	1 - 1.500 ml
Speed range	3.500 - 24.000 min ⁻¹
Motor rating input / output	500 / 300 W
Speed adjustment	infinitely variable
Speed display	scale
Dimensions (W x D x H)	65 x 80 x 240 mm
Weight	1,6 kg
Protection class acc. to DIN EN 60529	IP 20

DI 18 basic Ident. No. 3567000 472,00 EUR

Optional accessories:

Dispersing element S 18 N-19 G Ident. No. L004640 481,00 EUR



Package price
(DI 18 basic, S 18 N-10G, R 1826, R 182)
1.102,00 EUR
Package Ident. No. 9007900

Dispersing element S 18 N-10 G

With PTFE bearing

- Suitable for solvents and aggressive substances
- Dispersion tool can be sterilized using all methods




Material stainless steel	PTFE / AISI 316L
Diameter of stator	10 mm
Working range	1 - 100 ml
Max. temperature	180 °C

Ident. No. L004639 481,00 EUR

Plate stand R 1826

Stable, slip-free plate stand



Diameter of support rod	16 mm
Height	800 mm
(including PUR cover)	

Ident. No. 3160100 187,00 EUR

Boss head clamp R 182

Matches R 1826



Clamping range	
Stand	6 - 16 mm
Extension arm	6 - 16 mm

Ident. No. 2657700 40,00 EUR

DI 25 basic

Dispersing device

- High performance dispersion device for use in numerous applications
- Wide volume range
- Long service life

Volume range (H ₂ O)	1 - 2.000 ml
Speed range	8.000 - 24.000 min ⁻¹
Motor rating input / output	600 / 350 W
Speed adjustment	infinitely variable
Speed stability	1 %
Speed display	colour scale
Dimensions (W x D x H)	77 x 66 x 221 mm
Weight	1,6 kg
Protection class acc. to DIN EN 60529	IP 20

DI 25 basic Ident. No. 3130000 680,00 EUR

Optional accessories (further accessories on request):

Dispersing element S 25 N-10 G	Ident. No. 0594000	463,00 EUR
Dispersing element S 25 N-25 G	Ident. No. 1713300	644,00 EUR
Dispersing element S 25 N-25 F	Ident. No. 1713800	677,00 EUR



Package price
(DI 25 basic, S 25 N - 18 G, R 1826, R 182)
1.318,00 EUR
Package Ident. No. 9006400

Dispersing element S 25 N-18 G

With PTFE bearing

- Suitable for solvents and aggressive substances
- Dispersion tool can be sterilized using all methods




Material stainless steel	PTFE / AISI 316L
Diameter of stator	18 mm
Working range	10 - 1.500 ml
Max. temperature	180 °C

Ident. No. 0593400 502,00 EUR

Plate stand R 1826

Stable, slip-free plate stand



Diameter of support rod	16 mm
Height	800 mm
(including PUR cover)	

Ident. No. 3160100 187,00 EUR

Boss head clamp R 182

Matches R 1826



Clamping range	
Stand	6 - 16 mm
Extension arm	6 - 16 mm

Ident. No. 2657700 40,00 EUR